



**AP<sub>50/30</sub>**

## Real Time Intelligent Insufflator

- Utilizes real-time intracavitary pressure and flow sensing technology to deliver 100% CO<sub>2</sub> throughout the entire surgery
- Maintains constant pressure reading regardless of leaks or spikes
- One-way insufflator ensures no cross contamination of particles
- Eliminates bellowing effect

**Insuflow Port**

## Gas Conditioning Trocar

- Reduces lens fogging
- Improves visibility
- Protects patients from harmful effects of standard dry/cold or heated/dry Carbon Dioxide
- Dual lumen design provides uninterrupted flow regardless of instrument of scope present
- Available in 5mm, 8mm, 10mm, 12mm and bariatric sizes

# The LEXION System



**Insuflow**

## Heater Humidifier

- The original gas conditioning technology that heats and humidifies Carbon Dioxide to 95 °F/95% RH
- Provides clinical benefits by preventing intra-operative cooling and cell desiccation
- Insuflow® is designed to offer higher flow rates than standard tubing sets however for cases with heavy suction/smoke evacuation InsuflowPORT® should be used



**PneuVIEW XE**

## Continuous Smoke Evacuator

- .01 ULPA filtration removes hazardous smoke, gases, and harmful contaminants associated with combustion, including viral particles
- VFE/BFE verified
- Closed smoke evacuation continuously removed contaminants safely out of OR
- Custom flow settings with turbo button allows complete control of smoke evacuation at all times

**LEXION**

545 Atwater Circle, St. Paul, MN 55103  
877-953-9466 or 855-688-FLOW  
www.lexionmedical.com

Contact a Sales Representative today  
for a no-cost evaluation

AP 50/30 PIN 4100  
InsuflowPort® PIN 4100-5  
PneuVIEW XE® PIN 8660-XE  
©2019 LEXION Medical