


# URETERAL ACCESS SHEATH

With atraumatic tip-transition

## DESCRIPTION



### Ureteral Access Sheath

- ◆ Unique tip expanding system for atraumatic sheath transition
- ◆ Dilator lock system for secure connection between sheath and dilator
- ◆ Visual control of dilator fixation
- ◆ PTFE inner layer which provides smooth device passage
- ◆ Radiopaque sheath and dilator tip for easy positioning
- ◆ Hydrophilic coating to facilitate the introduction
- ◆ Stainless-steel coil reinforcement sheath which provides high anti-kink characteristics

Tube Length (cm)	Inner Diameter (CH/FR)	Outer Diameter (CH/FR)	Order Code
35	10	12	AS-1235
45	10	12	AS-1245
35	12	14	AS-1435
45	12	14	AS-1445
35	14	16	AS-1635
45	14	16	AS-1645