

Silicone Pessary Product Line

Bioteque offers a selection of soft, pliable, resilient pessaries in a variety of shapes and sizes, including malleable types, to relieve vaginal support abnormalities. Bioteque's patented pessaries are uniquely designed to easily fold for comfortable insertion and removal. All Bioteque pessaries include detailed physician instructions as well as patient information for self-care.



Ring

The Ring pessary, available with or without a supporting membrane, is primarily used for a mild first-degree uterine prolapse or a cystocele.

SKU(Without Support)	R#0	R#1	R#2	R#3	R#4	R#5	R#6	R#7	R#8	R#9
Dimension in inch	R1.75"	R2.00"	R2.25"	R2.50"	R2.75"	R3.00"	R3.25"	R3.50"	R3.75"	R4.00"
SKU(With Support)	Rs#0	Rs#1	Rs#2	Rs#3	Rs#4	Rs#5	Rs#6	Rs#7	Rs#8	Rs#9
Dimension in inch	R1.75S	R2.00S	R2.25S	R2.50S	R2.75S	R3.00S	R3.25S	R3.50S	R3.75S	R4.00S

(See back page for Ring with Support Fitting Set)



Cube

The Cube pessary, available with or without drains, may be used to relieve the symptoms of a cystocele or rectocele as well as a second or third-degree uterine prolapse.

SKU(Without Drains)	CU#0	CU#1	CU#2	CU#3	CU#4	CU#5	CU#6	CU#7	CU#8	CU#9	CU#10
Dimension in mm	CU25	CU29	CU33	CU37	CU41	CU45	CU50	CU56	CU63	CU70	CU75
SKU(With Drains)	CUd#0	CUd#1	CUd#2	CUd#3	CUd#4	CUd#5	CUd#6	CUd#7	CUd#8	CUd#9	CUd#10
Dimension in mm	CU25D	CU29D	CU33D	CU37D	CU41D	CU45D	CU50D	CU56D	CU63D	CU70D	CU75D

(See back page for Cube Fitting Set)



Gellhorn

The Gellhorn pessary with drains is used to relieve the symptoms of a cystocele or rectocele as well as a second or third-degree uterine prolapse. This pessary's stem with knob is easily folded to facilitate insertion and removal.

SKU(With Drains)	Gd#0	Gd#1	Gd#2	Gd#3	Gd#4	Gd#5	Gd#6	Gd#7	Gd#8
Dimension in inch	G1.50D	G1.75D	G2.00D	G2.25D	G2.50D	G2.75D	G3.00D	G3.25D	G3.50D

(See back page for Gellhorn Fitting Set)



Short-stem Gellhorn

The short-stem Gellhorn pessary with drains is also used to relieve the symptoms of a cystocele or rectocele as well as a second or third-degree uterine prolapse. It provides a better fit for some women.

SKU(With Drains)	GS#0	GS#1	GS#2	GS#3	GS#4	GS#5	GS#6	GS#7	GS#8
Dimension in inch	G1.50S	G1.75S	G2.00S	G2.25S	G2.50S	G2.75S	G3.00S	G3.25S	G3.50S



Donut

The Donut pessary may be used to relieve the stress of a second or third-degree uterine prolapse, a cystocele, or rectocele. This pessary, often used by post-menopausal women, is easily compressed for insertion and removal.

SKU	D#0	D#1	D#2	D#3	D#4	D#5	D#6	D#7
Dimension in inch	D2.00"	D2.25"	D2.50"	D2.75"	D3.00"	D3.25"	D3.50"	D3.75"



Dish

The Dish pessary, available with or without a supporting membrane, is used to relieve the stress of a second or third-degree uterine prolapse, a cystocele, or rectocele. This pessary may also relieve stress urinary incontinence.

SKU(Without Support)	DSH#0	DSH#1	DSH#2	DSH#3	DSH#4	DSH#5	DSH#6	DSH#7
Dimension in mm	DSH50	DSH55	DSH60	DSH65	DSH70	DSH75	DSH80	DSH85
SKU(With Support)	DSHs#0	DSHs#1	DSHs#2	DSHs#3	DSHs#4	DSHs#5	DSHs#6	DSHs#7
Dimension in mm	DSH50S	DSH55S	DSH60S	DSH65S	DSH70S	DSH75S	DSH80S	DSH85S

(See back page for Dish Fitting Set)



Oval (An unique product exclusive to Bioteque)

The Oval pessary, available with or without a supporting membrane, is commonly used for a mild first-degree uterine prolapse or a cystocele. This pessary's shape inhibits its rotation in women with a history of vaginal surgery.

SKU(Without Support)	OV#1	OV#2	OV#3	OV#4	OV#5	OV#6	OV#7	OV#8	OV#9
Dimension in inch	OV2.00"	OV2.25"	OV2.50"	OV2.75"	OV3.00"	OV3.25"	OV3.50"	OV3.75"	OV4.00"
SKU(With Support)	OVs#1	OVs#2	OVs#3	OVs#4	OVs#5	OVs#6	OVs#7	OVs#8	OVs#9
Dimension in inch	OV2.00S	OV2.25S	OV2.50S	OV2.75S	OV3.00S	OV3.25S	OV3.50S	OV3.75S	OV4.00S

(See back page for Oval Fitting Set)

Bioteque Pessary: CPT Code: A4562 Procedure Code: 57160

1-800-889-9008 Fax: (408) 526-9399

www.bioteque.com

Silicone Pessary Product Line



Marland (An unique product exclusive to Bioteque)

The Marland pessary, available with or without a supporting membrane, is used to relieve the symptoms of a second or third-degree uterine prolapse, a cystocele, or a rectocele. This pessary's elevated platform may also relieve stress urinary incontinence.

SKU(Without Support)	M#2	M#3	M#4	M#5	M#6	M#7	M#8
Dimension in inch	M2.25"	M2.50"	M2.75"	M3.00"	M3.25"	M3.50"	M3.75"
SKU(With Support)	Ms#2	Ms#3	Ms#4	Ms#5	Ms#6	Ms#7	Ms#8
Dimension in inch	M2.25S	M2.50S	M2.75S	M3.00S	M3.25S	M3.50S	M3.75S



Gehring

The malleable Gehring pessary with a support membrane is primarily used for a significant cystocele or rectocele and may also provide relief for a second or third-degree uterine prolapse. This pessary can be manually shaped to accommodate individual problems and anatomies.

SKU	GH#0	GH#1	GH#2	GH#3	GH#4	GH#5	GH#6	GH#7	GH#8	GH#9
Dimension in mm	GH50S	GH55S	GH60S	GH65S	GH70S	GH75S	GH80S	GH85S	GH90S	GH95S



Shatz

The Shatz pessary may relieve the symptoms of a cystocele as well as a second or third-degree uterine prolapse. Its convex shape ensures a snug fit while the drainage holes permit easy removal.

SKU	SH#0	SH#1	SH#2	SH#3	SH#4	SH#5	SH#6	SH#7	SH#8
Dimension in inch	SH1.50"	SH1.75"	SH2.00"	SH2.25"	SH2.50"	SH2.75"	SH3.00"	SH3.25"	SH3.50"



Hodge

The malleable Hodge pessary, available with or without a support membrane, is commonly used for a first or mild second-degree uterine prolapse and may also relieve a cystocele. These pessaries, designed to be manually shapeable for individual fitting, are usually prescribed for reproductive-aged women who prefer conservative therapy.

SKU(Without Support)	HD#0	HD#1	HD#2	HD#3	HD#4	HD#5	HD#6	HD#7	HD#8	HD#9
Dimension in mm	HD65	HD70	HD75	HD80	HD85	HD90	HD95	HD100	HD105	HD110
SKU(With Support)	HDs#0	HDs#1	HD#2	HD#3	HDs#4	HD#5	HDs#6	HD#7	HD#8	HDs#9
Dimension in mm	HD65S	HD70S	HD75S	HD80S	HD85S	HD90S	HD95S	HD100S	HD105	HD110S



Ring with Knob

The knobbed Ring pessary, available with or without a support membrane, has the same indications as the Ring and may also relieve stress urinary incontinence.

SKU(Without Support)	RK#0	RK#1	RK#2	RK#3	RK#4	RK#5	RK#6	RK#7
Dimension in inch	RK1.75"	RK2.00"	RK2.25"	RK2.50"	RK2.75"	RK3.00"	RK3.25"	RK3.50"
SKU(With Support)	RKs#0	RKs#1	RKs#2	RKs#3	RKs#4	RKs#5	RKs#6	RKs#7
Dimension in inch	RK1.75S	RK2.00S	RK2.25S	RK2.50S	RK2.75S	RK3.00S	RK3.25S	RK3.50S



Cup

The Cup pessary, available with or without a support membrane, is used for a first or a mild second-degree uterine prolapse.

SKU(Without Support)	CP#0	CP#1	CP#2	CP#3	CP#4	CP#5	CP#6	CP#7
Dimension in mm	CP50/18	CP55/23	CP60/28	CP65/33	CP70/38	CP75/43	CP80/48	CP85/52
SKU(With Support)	CPs#0	CPs#1	CPs#2	CPs#3	CPs#4	CPs#5	CPs#6	CPs#7
Dimension in mm	CP50/18S	CP55/23S	CP60/28S	CP65/33S	CP70/38S	CP75/43S	CP80/48S	CP85/52S



Silicone Inflatable Donut Pessary (An unique product exclusive to Bioteque)

This pessary can be used to relieve the symptoms of a mild second or third degree cystocele or a rectocele in combination with proidentia/prolapse.

SKU#	Dimension
ID 2	2.00"
ID 3	2.25"
ID 4	2.50"
ID 5	2.75"

Bioteque Pessary: CPT Code: A4562 Procedure Code: 57160

1-800-889-9008 Fax: (408) 526-9399

www.bioteque.com