


MAGNETIC BLACK-STAR Double-J Ureteral Stent-Set and Retrieval Device with magnet

SET COMPONENTS



Ureteral Stent

- ◆ Polyurethane (PUR)
- ◆ Drainage openings in pigtails and shaft
- ◆ Open tip
- ◆ Markings for exact positioning
- ◆ Radiopaque
- ◆ Magnet at bladder side



Pusher

- ◆ Length: 40 cm

Guide Wire

- ◆ PTFE-coated

Retrieval Device with magnet

- ◆ Soft PUR
- ◆ Easy tip orientation via marking line
- ◆ Tiemann tip, 30°
- ◆ Single packaging

URETERAL STENTS

MAGNETIC BLACK-STAR Double-J Ureteral Stent-Set and Retrieval Device with magnet

Open tip								
Size (FR/CH)	Length (cm)	Size Magnet of the Stent (FR/CH)	Order Code without Guide Wire with Retrieval Device	Order Code with Guide Wire with Retrieval Device	Guide Wire Length (cm)	Guide Wire Diameter	Guide Wire Tip	Guide Wire Stiffness
4.8	20	07	MR-414820	MR-364820	150	0.028"	straight	standard
4.8	22	07	MR-414822	MR-364822	150	0.028"	straight	standard
4.8	24	07	MR-414824	MR-364824	150	0.028"	straight	standard
4.8	26	07	MR-414826	MR-364826	150	0.028"	straight	standard
06	15	09	MR-410615	MR-360615	150	0.035"	straight	standard
06	18	09	MR-410618	MR-360618	150	0.035"	straight	standard
06	20	09	MR-414820	MR-364820	150	0.035"	straight	standard
06	22	09	MR-410622	MR-360622	150	0.035"	straight	standard
06	24	09	MR-410624	MR-360624	150	0.035"	straight	standard
06	26	09	MR-410626	MR-360626	150	0.035"	straight	standard
07	15	09	MR-410715	MR-360715	150	0.035"	straight	standard
07	18	09	MR-410718	MR-360718	150	0.035"	straight	standard
07	20	09	MR-410720	MR-360720	150	0.035"	straight	standard
07	22	09	MR-410722	MR-360722	150	0.035"	straight	standard
07	24	09	MR-410724	MR-360724	150	0.035"	straight	standard
07	26	09	MR-410726	MR-360726	150	0.035"	straight	standard

Offene Spitze ohne Retrieval Device								
Size (FR/CH)		Size Magnet of the Stent (FR/CH)	Order Code without Guide Wire without Retrieval Device	Order Code with Guide Wire without Retrieval Device	Guide Wire Length (cm)	Guide Wire Diameter	Guide Wire Tip	Guide Wire Stiffness
4.8	20	07	MR-404820	MR-304820	150	0.028"	straight	standard
4.8	22	07	MR-404822	MR-304822	150	0.028"	straight	standard
4.8	24	07	MR-404824	MR-304824	150	0.028"	straight	standard
4.8	26	07	MR-404826	MR-304826	150	0.028"	straight	standard
06	15	09	MR-400615	MR-300615	150	0.035"	straight	standard
06	18	09	MR-400618	MR-300618	150	0.035"	straight	standard
06	20	09	MR-400620	MR-300620	150	0.035"	straight	standard
06	22	09	MR-400622	MR-300622	150	0.035"	straight	standard
06	24	09	MR-400624	MR-300624	150	0.035"	straight	standard
06	26	09	MR-400626	MR-300626	150	0.035"	straight	standard
07	15	09	MR-400715	MR-300715	150	0.035"	straight	standard
07	18	09	MR-400718	MR-300718	150	0.035"	straight	standard
07	20	09	MR-400720	MR-300720	150	0.035"	straight	standard
07	22	09	MR-400722	MR-300722	150	0.035"	straight	standard
07	24	09	MR-400724	MR-300724	150	0.035"	straight	standard
07	26	09	MR-400726	MR-300726	150	0.035"	straight	standard

Retrieval Device

Size (FR/CH)	Length (cm)	Order Code suitable for FR/CH 4.8	Order Code suitable for FR/CH 6/7
09	40	R-Device-S	-
15	40	-	R-Device-L