

URETERAL STENTS


GREEN-STAR



GREEN-STAR Double-J Ureteral Stent-Set

Soft material

SET COMPONENTS



Ureteral Stent

- ♦ Polyurethane (PUR)
- ♦ Special soft material for high patient comfort
- ♦ Drainage openings in pigtails and shaft
- ♦ Open or closed tip
- ♦ Markings for exact positioning
- ♦ High radiopaque visibility

Pusher

- ♦ Length: 40 cm

Guide Wire

- ♦ PTFE-coated

Clamp

- ♦ Only in sets with closed tip



GREEN-STAR Double-J Ureteral Stent-Set

Soft material


Open tip							
Size (FR/CH)	Length (cm)	Order Code without Guide Wire	Order Code with Guide Wire	Guide Wire Length (cm)	Guide Wire Diameter	Guide Wire Tip	Guide Wire Stiffness
06	24	EG-410624	EG-360624	125	0.035"	straight	standard
06	26	EG-410626	EG-360626	125	0.035"	straight	standard
06	28	EG-410628	EG-360628	125	0.035"	straight	standard
06	30	EG-410630	EG-360630	125	0.035"	straight	standard
07	24	EG-410724	EG-360724	125	0.035"	straight	standard
07	26	EG-410726	EG-360726	125	0.035"	straight	standard
07	28	EG-410728	EG-360728	125	0.035"	straight	standard
07	30	EG-410730	EG-360730	125	0.035"	straight	standard
08	24		EG-360824	125	0.035"	straight	standard
08	26	EG-410826	EG-360826	125	0.035"	straight	standard
08	28	EG-410828	EG-360828	125	0.035"	straight	standard
08	30		EG-360830	125	0.035"	straight	standard

Closed tip							
Size (FR/CH)	Length (cm)	Order Code without Guide Wire	Order Code with Guide Wire	Guide Wire Length (cm)	Guide Wire Diameter	Guide Wire Tip	Guide Wire Stiffness
06	24	EG-400624	EG-390624	125	0.035"	straight	standard
06	26	EG-400626	EG-390626	125	0.035"	straight	standard
06	28	EG-400628	EG-390628	125	0.035"	straight	standard
06	30	EG-400630	EG-390630	125	0.035"	straight	standard
07	24		EG-390724	125	0.035"	straight	standard
07	26	EG-400726	EG-390726	125	0.035"	straight	standard
07	28	EG-400728	EG-390728	125	0.035"	straight	standard
07	30	EG-400730	EG-390730	125	0.035"	straight	standard
08	26	EG-400826	EG-390826	125	0.035"	straight	standard
08	28	EG-400828	EG-390828	125	0.035"	straight	standard

GREEN-STAR Double-J Ureteral Stent-Set


Soft material, flexible length

SET COMPONENTS



Ureteral Stent

- ♦ Polyurethane (PUR)
- ♦ Special soft material for high patient comfort
- ♦ Multiple coils on both ends for flexible shaft length from 22 to 30 cm
- ♦ Drainage openings in pigtails and shaft
- ♦ Open or closed tip
- ♦ Markings for exact positioning
- ♦ High radiopaque visibility



Pusher

- ♦ Length: 40 cm

Guide Wire

- ♦ PTFE-coated

Clamp

- ♦ Only in sets with closed tip



GREEN-STAR Double-J Ureteral Stent-Set

Soft material, flexible length


Open tip							
Size (FR/CH)	Length (cm)	Order Code without Guide Wire	Order Code with Guide Wire	Guide Wire Length (cm)	Guide Wire Diameter	Guide Wire Tip	Guide Wire Stiffness
06	22-30	EG-410622-30-BS	EG-360622-30-BS	125	0.035"	straight	standard
07	22-30	EG-410722-30-BS	EG-360722-30-BS	125	0.035"	straight	standard
08	22-30	EG-410822-30-BS	EG-360822-30-BS	125	0.035"	straight	standard

Closed tip							
Size (FR/CH)	Length (cm)		Order Code with Guide Wire	Guide Wire Length (cm)	Guide Wire Diameter	Guide Wire Tip	Guide Wire Stiffness
06	22-30		EG-390622-30-BS	125	0.035"	straight	standard
07	22-30		EG-390722-30-BS	125	0.035"	straight	standard

GREEN-STAR Double-J Ureteral Stent-Set

Soft material, steerable

SET COMPONENTS



Ureteral Stent

- ♦ Polyurethane (PUR)
- ♦ Special soft material for high patient comfort
- ♦ Drainage openings in pigtails and shaft
- ♦ Open or closed tip
- ♦ Markings for exact positioning
- ♦ High radiopaque visibility

Steering Device

- ♦ Pusher
- ♦ Disconnecter
- ♦ Length: 53 cm

Guide Wire

- ♦ PTFE-coated

Clamp

- ♦ Only in sets with closed tip



GREEN-STAR Double-J Ureteral Stent-Set Soft material, steerable


Open tip							
Size (FR/CH)	Length (cm)	Order Code without Guide Wire	Order Code with Guide Wire	Guide Wire Length (cm)	Guide Wire Diameter	Guide Wire Tip	Guide Wire Stiffness
06	24	EG-470624	EG-480624	150	0.035"	straight	standard / movable core
06	26	EG-470626	EG-480626	150	0.035"	straight	standard / movable core
06	28	EG-470628	EG-480628	150	0.035"	straight	standard / movable core
06	30	EG-470630	EG-480630	150	0.035"	straight	standard / movable core
07	26	EG-470726	EG-480726	150	0.035"	straight	standard / moveable core
07	28	EG-470728	EG-480728	150	0.035"	straight	standard / moveable core
07	30	EG-470730	EG-480730	150	0.035"	straight	standard / moveable core
08	24		EG-480824	150	0.035"	straight	standard / moveable core
08	26	EG-470826		150	0.035"	straight	standard / moveable core
08	28	EG-470828		150	0.035"	straight	standard / moveable core
08	30	EG-470830	EG-480830	150	0.035"	straight	standard / moveable core

Closed tip							
Size (FR/CH)	Length (cm)		Order Code with Guide Wire	Guide Wire Length (cm)	Guide Wire Diameter	Guide Wire Tip	Guide Wire Stiffness
06	24		EG-490624	125	0.035"	J-tip	standard / movable core
06	26		EG-490626	125	0.035"	J-tip	standard / movable core
06	28		EG-490628	125	0.035"	J-tip	standard / movable core
07	28		EG-490728	125	0.035"	J-tip	standard / movable core
08	26		EG-490826	125	0.035"	J-Tip	standard / movable core
08	28		EG-490828	125	0.035"	J-tip	standard / movable core

GREEN-STAR Double-J Ureteral Stent-Set

Soft material, steerable with coupling

SET COMPONENTS



Ureteral Stent

- ♦ Polyurethane (PUR)
- ♦ Special soft material for high patient comfort
- ♦ Drainage openings in pigtails and shaft
- ♦ Open or closed tip
- ♦ Markings for exact positioning
- ♦ High radiopaque visibility

Steering Device with Coupling

- ♦ Pusher
- ♦ Disconnecter
- ♦ Coupling to fasten guide wire
- ♦ Length: 43 cm

Guide Wire

- ♦ PTFE-coated

GREEN-STAR Double-J Ureteral Stent-Set Soft material, steerable with coupling


Open tip							
Size (FR/CH)	Length (cm)	Order Code without Guide Wire	Order Code with Guide Wire	Guide Wire Length (cm)	Guide Wire Diameter	Guide Wire Tip	Guide Wire Stiffness
06	24	EG-570624	EG-580624	150	0.035"	straight	standard / movable core
06	26	EG-570626	EG-580626	150	0.035"	straight	standard / movable core
06	28	EG-570628	EG-580628	150	0.035"	straight	standard / movable core
06	30		EG-580630	150	0.035"	straight	standard / movable core
07	24	EG-570724	EG-580724	150	0.035"	straight	standard / movable core
07	26	EG-570726	EG-580726	150	0.035"	straight	standard / movable core
07	28	EG-570728	EG-580728	150	0.035"	straight	standard / movable core
07	30	EG-570730	EG-580730	150	0.035"	straight	standard / movable core

Closed tip							
Size (FR/CH)	Length (cm)		Order Code with Guide Wire	Guide Wire Length (cm)	Guide Wire Diameter	Guide Wire Tip	Guide Wire Stiffness
06	24		EG-590624	125	0.035"	J-tip	standard / movable core
06	26		EG-590626	125	0.035"	J-tip	standard / movable core
06	28		EG-590628	125	0.035"	J-tip	standard / movable core
06	30		EG-590630	125	0.035"	J-tip	standard / movable core
07	26		EG-590726	125	0.035"	J-tip	standard / movable core
07	28		EG-590728	125	0.035"	J-tip	standard / movable core
08	24		EG-590824	125	0.035"	J-tip	standard / movable core
08	26		EG-590826	125	0.035"	J-tip	standard / movable core
08	28		EG-590828	125	0.035"	J-tip	standard / movable core

GREEN-STAR MONO Mono-J Ureteral Stent-Set

Soft material

SET COMPONENTS



Ureteral Stent

- ♦ Polyurethane (PUR)
- ♦ Special soft material for high patient comfort
- ♦ Drainage openings only in pigtail
- ♦ Open or closed tip
- ♦ Markings for exact positioning
- ♦ High radiopaque visibility

Guide Wire

- ♦ PTFE-coated

Adapter

- ♦ For attachment of the urinary drainage bag

Clamp

- ♦ Only in sets with closed tip



GREEN-STAR MONO Mono-J Ureteral Stent-Set

Soft material

Open tip						
Size (FR/CH)	Length (cm)	Order Code with Guide Wire	Guide Wire Length (cm)	Guide Wire Diameter	Guide Wire Tip	Guide Wire Stiffness
06	70	EG-350671	175	0.035"	straight	standard
07	70	EG-350771	175	0.035"	straight	standard
08	70	EG-350871	175	0.038"	straight	standard

Closed tip						
Size (FR/CH)	Length (cm)	Order Code with Guide Wire	Guide Wire Length (cm)	Guide Wire Diameter	Guide Wire Tip	Guide Wire Stiffness
06	70	EG-350670	125	0.035"	straight	standard
07	70	EG-350770	125	0.035"	straight	standard
08	70	EG-350870	125	0.035"	straight	standard