



The Lipo Cart is a complete and powerful Liposuction cart for tumescent assisted liposuction.

The Lipo Cart unites all high-quality equipment with simple and easy operation at highest flexibility.

Power Liposuction

Infiltration

Vacuum

Lipo Cart

Complete Liposuction system

The Lipo Cart represents our fine solution as a complete power assisted system for Liposuction.

Our product portfolio for liposuction features a power assisted device that is easier for surgeons to control, allowing more precise and thorough fat removal. Nouvag's LipoSurg produces rapid, small reciprocating vibrations to break up fat, so that it can easily be sucked out through the cannula.

Surgeons can remove more fat in a shorter period of time, which leads to enhanced results, a safer procedure, and faster recovery. Liposuction is performed using the tumescent liposuction technique, in which a liquid solution is delivered to create an effect of fat tunneling, and so enables easier fat transfer through the cannula and out of the body. The infiltration of an anesthetic substance into the treatment area is done by the LipoSurg integrated sophisticated infiltration system.