

LipoSurg

Vibration liposuction with tumescent infiltration

The LipoSurg system combines the functions of infiltration and motor-assisted vibration liposuction into a single, compact device. The desired values for infiltration and liposuction are set by selecting the function and pressing the ascending or descending arrow keys. Once the values are set, the liposuction or infiltration can be controlled by the VARIO or the ON/OFF pedal.

The extraction of fat requires the connection of this compact system with the Nouvag high vacuum pump "Vacuson 60". The high flow rate of 60 liters per minute is ideal for efficient liposuction.

Integrated high volume infiltration system with balanced pressure for tissue protection and a flow rate of 400 ml/min.

Lipo Surg



400 ml/min infiltration

up to 4200 strokes

Power Liposuction