

# Liposuction cannulas

High quality finished, cell friendly soft tissue cannulas from non-aggressive sculpting and extraction up to aggressive cannulas for large extraction volume. Especially designed cannulas for the use of tumescent local anesthesia power liposuction technique.



**Straight Entner cannulas**  
for abdominal surgery

| Cannula-Ø, mm | Length mm | No. of holes | Hole-Ø mm | REF  |
|---------------|-----------|--------------|-----------|------|
| 1.5           | 150       | 1            | oval      | 4361 |
| 2.0           | 150       | 1            | oval      | 4364 |
| 2.0           | 150       | 18           | 1.0       | 4373 |
| 3.0           | 150       | 18           | 1.5       | 4374 |
| 3.0           | 200       | 22           | 1.5       | 4378 |
| 3.0           | 300       | 30           | 1.5       | 4387 |
| 4.0           | 200       | 22           | 2.0       | 4379 |
| 4.0           | 300       | 30           | 1.5       | 4388 |

**Curved cannulas**  
for facial and neck liposuction

|     |     |    |     |      |
|-----|-----|----|-----|------|
| 3.0 | 200 | 22 | 1.5 | 4381 |
|-----|-----|----|-----|------|

**Bayonet cannulas**  
for thigh and lower leg liposuction

|     |     |    |     |      |
|-----|-----|----|-----|------|
| 3.0 | 200 | 22 | 1.5 | 4362 |
| 3.0 | 300 | 30 | 1.5 | 4365 |
| 4.0 | 200 | 22 | 1.5 | 4368 |
| 4.0 | 300 | 30 | 1.5 | 4372 |

**Mercedes**  
most popular for extraction



| Cannula-Ø mm | Length mm | REF  |
|--------------|-----------|------|
| 3.0          | 150       | 4487 |
| 3.0          | 250       | 4488 |
| 4.0          | 150       | 4489 |
| 4.0          | 250       | 4490 |

**Toledo**



|     |     |      |
|-----|-----|------|
| 2.0 | 50  | 4595 |
| 3.0 | 150 | 4596 |
| 4.0 | 250 | 4597 |
| 5.0 | 300 | 4598 |

**One central, two lateral**



| Cannula-Ø mm | Length mm | REF  |
|--------------|-----------|------|
| 3.0          | 200       | 4480 |
| 3.0          | 300       | 4481 |
| 4.0          | 300       | 4482 |
| 5.0          | 300       | 4483 |
| 6.0          | 300       | 4484 |

**One lateral**



|     |     |      |
|-----|-----|------|
| 2.0 | 150 | 4364 |
| 3.0 | 200 | 4485 |
| 3.0 | 300 | 4486 |



## Autologous fat transfer

An autologous fat transfer has many applications in plastic and aesthetic surgery. The autologous fat cells of the patient are harvested by gentle liposuction.

It is crucial that no damage occurs to the delicate fat cells during the

harvesting and transplanting process. Therefore highly quality precision Instrumentation in combination with a dedicated powered liposuction unit is absolutely vital.



**REF 4491:** Liposuction cannula for the harvesting of autologous fat cells. 15 cm length, Ø 2,5 mm, 4 perforations, with Luer-Lock connection.



**REF 4498:** Liposuction cannula for harvesting of autologous fat cells. 18 cm length, Ø 3 mm, sharp perforations, with Luer-Lock connection.