


RENODRAIN-INTEGRALE Nephrostomy Puncture- and Exchange Set

Seldinger technique, 100% silicone, integrated balloon

SET COMPONENTS



Balloon Catheter

- ◆ Silicone for optimal tissue compatibility
- ◆ Opening behind balloon to provide calyce drainage
- ◆ Integrated balloon
- ◆ Radiopaque contrast stripes
- ◆ Ultrasonic visible reinforced tip
- ◆ Markings for exact positioning



Puncture Needle

- ◆ 2-part needle
- ◆ Medical stainless steel
- ◆ 18 G, length: 20 cm
- ◆ Ultrasonic visible

Guide Wire (Lunderquist)

- ◆ PTFE-coated

RENODRAIN-INTEGRALE Nephrostomy Puncture- and Exchange Set

Seldinger technique, 100% silicone, integrated balloon

Single Catheter			
Size (FR/CH)	Length (cm)	Balloon (ml)	Order Code Single Catheter
10	40	1.5	NE-261005
12	30	2.0	NE-261225
12	40	2.0	NE-261205
14	30	3.0	NE-261425
14	40	3.0	NE-261405
16	30	5.0	NE-261625
16	40	5.0	NE-261605
18	30	5.0	NE-261825
18	40	5.0	NE-261805
20	40	5.0	NE-262005
22	40	5.0	NE-262205

Puncture Set								
Size (FR/CH)	Length (cm)	Balloon (ml)	Order Code Puncture Set	Dilators (FR/CH) *with split sheath	Guide Wire Length (cm)	Guide Wire Diameter	Guide Wire Tip (r=1.5 mm)	Guide Wire Stiffness
10	40	1.5	RE-401033	FR/CH 10/12*	70	0.025"	J-Tip	stiff
12	40	2.0	RE-401233	FR/CH 10/14*	70	0.025"	J-Tip	stiff
14	40	3.0	RE-401433	FR/CH 10/14/16*	80	0.035"	J-Tip	stiff
16	40	5.0	RE-401633	FR/CH 12/14/16*	80	0.035"	J-Tip	stiff

Exchange Set								
Size (FR/CH)	Length (cm)	Balloon (ml)	Order Code Exchange Set	Dilators (FR/CH) *with split sheath	Guide Wire Length (cm)	Guide Wire Diameter	Guide Wire Tip	Guide Wire Stiffness
10	40	1.5	RE-401032	FR/CH 10/12*	80	0.035"	straight	stiff
12	40	2.0	RE-401232	FR/CH 12/14*	80	0.035"	straight	stiff
14	40	3.0	RE-401432	FR/CH 14/16*	80	0.035"	straight	stiff
16	40	5.0	RE-401632	FR/CH 16/18*	80	0.035"	straight	stiff