



RENODRAIN-YELLOW Nephrostomy Puncture- and Exchange Set Seldinger technique, Treated with PC, with pigtail

SET COMPONENTS



Nephrostomy Catheter

- ◆ Polyurethane (PUR)
- ◆ Treated with PC to minimise the risk of infection and encrustation, therefore extended indwelling time
- ◆ Drainage openings in pigtail
- ◆ Luer Lock Connector



Puncture Needle

- ◆ 2-part needle
- ◆ Medical stainless steel
- ◆ 18 G, length: 20 cm
- ◆ Ultrasonic visible

Guide Wire (Lunderquist)

- ◆ PTFE-coated

Universal Adapter

- ◆ Male Adapter with Luer Lock Connector
- ◆ 10 cm Silicone tube

NEPHROSTOMY

RENODRAIN-YELLOW Nephrostomy Puncture- and Exchange Set Seldinger technique, Treated with PC, with pigtail

Single Catheter		
Size (FR/CH)	Length (cm)	Order Code Pigtail Catheter
06	24	RP-410640
06	34	RP-410650
07	24	RP-410740
07	34	RP-410750
08	24	RP-410840
08	34	RP-410850
09	24	RP-410940
09	34	RP-410950
10.5	30	RP-411050
12	30	RP-411250
14	28	RP-411440

Puncture Set							
Size (FR/CH)	Length (cm)	Order Code Puncture Set	Dilators (FR/CH)	Guide Wire Length (cm)	Guide Wire Diameter	Guide Wire Tip (r=1.5 mm)	Guide Wire Stiffness
06	24	RP-400640	05/07	80	0.035"	J-Tip	stiff
06	34	RP-400650	05/07	80	0.035"	J-Tip	stiff
07	24	RP-400740	06/08	80	0.035"	J-Tip	stiff
07	34	RP-400750	06/08	80	0.035"	J-Tip	stiff
08	24	RP-400840	07/09	80	0.035"	J-Tip	stiff
08	34	RP-400850	07/09	80	0.035"	J-Tip	stiff
09	24	RP-400940	08/10	80	0.035"	J-Tip	stiff
09	34	RP-400950	08/10	80	0.035"	J-Tip	stiff
10.5	30	RP-401050	08/12	80	0.035"	J-Tip	stiff
12	30	RP-401250	10/14	80	0.035"	J-Tip	stiff
14	28	RP-401440	10/14/16	80	0.035"	J-Tip	stiff

Exchange Set							
Size (FR/CH)	Length (cm)	Order Code Exchange Set	Dilators (FR/CH)	Guide Wire Length (cm)	Guide Wire Diameter	Guide Wire Tip	Guide Wire Stiffness
06	24	RP-400642	06/07	80	0.035"	straight	stiff
06	34	RP-400652	06/07	80	0.035"	straight	stiff
07	24	RP-400742	07/08	80	0.035"	straight	stiff
07	34	RP-400752	07/08	80	0.035"	straight	stiff
08	24	RP-400842	08/09	80	0.035"	straight	stiff
08	34	RP-400852	08/09	80	0.035"	straight	stiff
09	24	RP-400942	09/10	80	0.035"	straight	stiff
09	34	RP-400952	09/10	80	0.035"	straight	stiff
10.5	30	RP-401052	10/12	80	0.035"	straight	stiff
12	30	RP-401252	12/14	80	0.035"	straight	stiff
14	28	RP-401442	12/14/16	80	0.035"	straight	stiff