

Vacuson 60

High performance vacuum pump, ideal for liposuction

The suction capacity of 60 liters per minute predestines this pump for the liposuction environment. Large volumes of fat can be removed in a short time.

Depending on the requirement, two 2 liter secretion bottles or two 5 liter secretion bottles can each be hung on the bottle holders of the Vacuson 60 or the VEXIO Cart device trolley. With the MONOKIT system, it is also possible to work with disposable pouches and thus reduce the cleaning and reprocessing effort.



REF 4237, Vacuson 60 with pedal

- 115 V version
- 230 V version

60 l/min