

VEXIO Cart "Profi"

Complete, safe, flexible



This robust and easy maneuverable device cart is the base for an uncomplicated, man and machine friendly accommodation of your device combination.

At the place of use, the VEXIO Cart is connected to the grid. As a result, all devices on the cart are ready for use immediately. The cables of the devices used disappear in the built in cable channel.

Integrated power supply with 3 sockets for maximum maneuverability

Antistatic, individually lockable rollers ensure high security

Cable ducts in the center column ensure order and contribute to safe handling

REF 1898: VEXIO Cart "Profi"

REF 1894: Bottle holder is optional