


URETERAL STENTS

WHITE-STAR



WHITE-STAR Double-J Ureteral Stent-Set

SET COMPONENTS



Ureteral Stent

- ♦ Polyurethane (PUR)
- ♦ Drainage openings in pigtails and shaft
- ♦ Open or closed tip
- ♦ Markings for exact positioning
- ♦ Radiopaque

Pusher

- ♦ Length: 40 cm

Guide Wire

- ♦ PTFE-coated

Clamp

- ♦ 2 clamps in sets with closed tip



WHITE-STAR Double-J Ureteral Stent-Set


Open tip							
Size (FR/CH)	Length (cm)	Order Code without Guide Wire	Order Code with Guide Wire	Guide Wire Length (cm)	Guide Wire Diameter	Guide Wire Tip	Guide Wire Stiffness
4.7	24	EN-410524	EN-360524	125	0.028"	straight	standard
4.7	26	EN-410526	EN-360526	125	0.028"	straight	standard
06	24	EN-410624	EN-360624	125	0.035"	straight	standard
06	26	EN-410626	EN-360626	125	0.035"	straight	standard
06	28	EN-410628	EN-360628	125	0.035"	straight	standard
06	30	EN-410630	EN-360630	125	0.035"	straight	standard
07	24	EN-410724	EN-360724	125	0.035"	straight	standard
07	26	EN-410726	EN-360726	125	0.035"	straight	standard
07	28	EN-410728	EN-360728	125	0.035"	straight	standard
07	30	EN-410730	EN-360730	125	0.035"	straight	standard
08	24	EN-410824	EN-360824	125	0.035"	straight	standard
08	26	EN-410826	EN-360826	125	0.035"	straight	standard
08	28	EN-410828	EN-360828	125	0.035"	straight	standard
08	30	EN-410830	EN-360830	125	0.035"	straight	standard

Closed tip							
Size (FR/CH)	Length (cm)	Order Code without Guide Wire	Order Code with Guide Wire	Guide Wire Length (cm)	Guide Wire Diameter	Guide Wire Tip	Guide Wire Stiffness
4.7	24	EN-400524	EN-300524	125	0.028"	straight	standard
4.7	26	EN-400526	EN-300526	125	0.028"	straight	standard
06	24	EN-400624	EN-300624	125	0.035"	straight	standard
06	26	EN-400626	EN-300626	125	0.035"	straight	standard
06	28	EN-400628	EN-300628	125	0.035"	straight	standard
06	30	EN-400630	EN-300630	125	0.035"	straight	standard
07	24	EN-400724	EN-300724	125	0.035"	straight	standard
07	26	EN-400726	EN-300726	125	0.035"	straight	standard
07	28	EN-400728	EN-300728	125	0.035"	straight	standard
07	30	EN-400730	EN-300730	125	0.035"	straight	standard
08	24	EN-400824	EN-300824	125	0.035"	straight	standard
08	26	EN-400826	EN-300826	125	0.035"	straight	standard
08	28	EN-400828	EN-300828	125	0.035"	straight	standard
08	30		EN-300830	125	0.035"	straight	standard

WHITE-STAR Double-J Ureteral Stent-Set


Flexible length

SET COMPONENTS



Ureteral Stent

- ♦ Polyurethane (PUR)
- ♦ Multiple coils on both ends for flexible shaft length from 22 to 30 cm
- ♦ Less risk of dislocation
- ♦ Drainage openings in pigtails and shaft
- ♦ Open or closed tip
- ♦ Markings for exact positioning
- ♦ Radiopaque



Pusher

- ♦ Length: 40 cm

Guide Wire

- ♦ PTFE-coated

Clamp

- ♦ 2 clamps in sets with closed tip



WHITE-STAR Double-J Ureteral Stent-Set


Flexible length

Open tip							
Size (FR/CH)	Length (cm)	Order Code without Guide Wire	Order Code with Guide Wire	Guide Wire Length (cm)	Guide Wire Diameter	Guide Wire Tip	Guide Wire Stiffness
06	22-30	EN-410622-30-BS	EN-360622-30-BS	125	0.035"	straight	standard
07	22-30	EN-410722-30-BS	EN-360722-30-BS	125	0.035"	straight	standard
08	22-30	EN-410822-30-BS	EN-360822-30-BS	125	0.035"	straight	standard

Closed tip							
Size (FR/CH)	Length (cm)		Order Code with Guide Wire	Guide Wire Length (cm)	Guide Wire Diameter	Guide Wire Tip	Guide Wire Stiffness
06	22-30		EN-390622-30-BS	125	0.035"	straight	standard
07	22-30		EN-390722-30-BS	125	0.035"	straight	standard
08	22-30		EN-390822-30-BS	125	0.035"	straight	standard

WHITE-STAR Double-J Ureteral Stent-Set Steerable

SET COMPONENTS



Ureteral Stent

- ♦ Polyurethane (PUR)
- ♦ Drainage openings in pigtails and shaft
- ♦ Open or closed tip
- ♦ Markings for exact positioning
- ♦ Radiopaque

Steering Device

- ♦ Pusher
- ♦ Disconnecter
- ♦ Length: 53 cm

Guide Wire

- ♦ PTFE-coated

Clamp

- ♦ 2 clamps in sets with closed tip



WHITE-STAR Double-J Ureteral Stent-Set

Steerable


Open tip							
Size (FR/CH)	Length (cm)	Order Code without Guide Wire	Order Code with Guide Wire	Guide Wire Length (cm)	Guide Wire Diameter	Guide Wire Tip	Guide Wire Stiffness
4.7	24	ES-360524	ES-370524	145	0.028"	straight	standard
4.7	26	ES-360526	ES-370526	145	0.028"	straight	standard
06	24	ES-360624	ES-370624	150	0.035"	straight	standard / movable core
06	26	ES-360626	ES-370626	150	0.035"	straight	standard / movable core
06	28	ES-360628	ES-370628	150	0.035"	straight	standard / movable core
06	30	ES-360630	ES-370630	150	0.035"	straight	standard / movable core
07	24	ES-360724	ES-370724	150	0.035"	straight	standard / movable core
07	26	ES-360726	ES-370726	150	0.035"	straight	standard / movable core
07	28	ES-360728	ES-370728	150	0.035"	straight	standard / movable core
07	30	ES-360730	ES-370730	150	0.035"	straight	standard / movable core
08	24		ES-370824	150	0.035"	straight	standard / movable core
08	26	ES-360826	ES-370826	150	0.035"	straight	standard / movable core
08	28	ES-360828	ES-370828	150	0.035"	straight	standard / movable core
08	30	ES-360830	ES-370830	150	0.035"	straight	standard / movable core

Closed tip							
Size (FR/CH)	Length (cm)		Order Code with Guide Wire	Guide Wire Length (cm)	Guide Wire Diameter	Guide Wire Tip	Guide Wire Stiffness
4.7	24		ES-380524	125	0.028"	straight	standard
4.7	26		ES-380526	125	0.028"	straight	standard
06	24		ES-380624	125	0.035"	J-tip	standard / movable core
06	26		ES-380626	125	0.035"	J-tip	standard / movable core
06	28		ES-380628	125	0.035"	J-tip	standard / movable core
06	30		ES-380630	125	0.035"	J-tip	standard / movable core
07	26		ES-380726	125	0.035"	J-tip	standard / movable core
07	28		ES-380728	125	0.035"	J-tip	standard / movable core
07	30		ES-380730	125	0.035"	J-tip	standard / movable core

WHITE-STAR Double-J Ureteral Stent-Set

Steerable with coupling

SET COMPONENTS



Ureteral Stent

- ♦ Polyurethane (PUR)
- ♦ Drainage openings in pigtails and shaft
- ♦ Open or closed tip
- ♦ Markings for exact positioning
- ♦ Radiopaque

Steering Device with Coupling

- ♦ Pusher
- ♦ Disconnecter
- ♦ Coupling to fasten guide wire
- ♦ Length: 43 cm

Guide Wire

- ♦ PTFE-coated

WHITE-STAR Double-J Ureteral Stent-Set


Steerable with coupling

Open tip							
Size (FR/CH)	Length (cm)	Order Code without Guide Wire	Order Code with Guide Wire	Guide Wire Length (cm)	Guide Wire Diameter	Guide Wire Tip	Guide Wire Stiffness
4.7	24	ES-570524	ES-580524	145	0.028"	straight	standard
4.7	26	ES-570526	ES-580526	145	0.028"	straight	standard
06	24		ES-580624	150	0.035"	straight	standard / movable core
06	26	ES-570626	ES-580626	150	0.035"	straight	standard / movable core
06	28	ES-570628	ES-580628	150	0.035"	straight	standard / movable core
06	30	ES-570630	ES-580630	150	0.035"	straight	standard / movable core
07	24	ES-570724		150	0.035"	straight	standard / movable core
07	26	ES-570726	ES-580726	150	0.035"	straight	standard / movable core
07	28	ES-570728	ES-580728	150	0.035"	straight	standard / movable core
07	30	ES-570730	ES-580730	150	0.035"	straight	standard / movable core
08	24		ES-580824	150	0.035"	straight	standard / movable core
08	26	ES-570826	ES-580826	150	0.035"	straight	standard / movable core
08	28		ES-580828	150	0.035"	straight	standard / movable core
08	30		ES-580830	150	0.035"	straight	standard / movable core

Closed tip							
Size (FR/CH)	Length (cm)		Order Code with Guide Wire	Guide Wire Length (cm)	Guide Wire Diameter	Guide Wire Tip	Guide Wire Stiffness
4.7	24		ES-590524	125	0.028"	straight	standard
4.7	26		ES-590526	125	0.028"	straight	standard
06	24		ES-590624	125	0.035"	J-tip	standard / movable core
06	26		ES-590626	125	0.035"	J-tip	standard / movable core
06	28		ES-590628	125	0.035"	J-tip	standard / movable core
06	30		ES-590630	125	0.035"	J-tip	standard / movable core
07	24		ES-590724	125	0.035"	J-tip	standard / movable core
07	26		ES-590726	125	0.035"	J-tip	standard / movable core
07	28		ES-590728	125	0.035"	J-tip	standard / movable core
07	30		ES-590730	125	0.035"	J-tip	standard / movable core
08	24		ES-590824	125	0.035"	J-tip	standard / movable core
08	26		ES-590826	125	0.035"	J-tip	standard / movable core

WHITE-STAR MONO Mono-J Ureteral Stent-Set

SET COMPONENTS



Ureteral Stent

- ♦ Polyurethane (PUR)
- ♦ Drainage openings only in pigtail
- ♦ Open or closed tip
- ♦ Markings for exact positioning
- ♦ Radiopaque

Guide Wire

- ♦ PTFE-coated

Adapter

- ♦ For attachment of the urinary drainage bag

Clamp

- ♦ Only in sets with closed tip



WHITE-STAR MONO Mono-J Ureteral Stent-Set

Open tip						
Size (FR/CH)	Length (cm)	Order Code with Guide Wire	Guide Wire Length (cm)	Guide Wire Diameter	Guide Wire Tip	Guide Wire Stiffness
04	70	EN-350471	125	0.025"	straight	standard
4.7	70	EN-350571	175	0.028"	straight	standard
06	70	EN-350671	175	0.035"	straight	standard
07	70	EN-350771	175	0.035"	straight	standard
08	70	EN-350871	175	0.038"	straight	standard
09	70	EN-350971	175	0.045"	straight	standard
10	70	EN-351071	175	0.045"	straight	standard
12	70	EN-351271	175	0.045"	straight	standard

Closed tip						
Size (FR/CH)	Length (cm)	Order Code with Guide Wire	Guide Wire Length (cm)	Guide Wire Diameter	Guide Wire Tip	Guide Wire Stiffness
04	70	EN-350470	125	0.025"	straight	standard
4.7	70	EN-350570	125	0.028"	straight	standard
06	70	EN-350670	125	0.035"	straight	standard
07	70	EN-350770	125	0.035"	straight	standard
08	70	EN-350870	125	0.035"	straight	standard
09	70	EN-350970	125	0.045"	straight	standard
10	70	EN-351070	125	0.045"	straight	standard
12	70	EN-351270	125	0.045"	straight	standard